



■ Consumables for CyBi®-FeliX Systems

Pipetting tips for CyBi®-FeliX systems are available in packaging units (PU) with pre-racked CyBi®-RoboTipTrays or CyBi®-TipBoxes.

CHOICE™ Pipetting Head

using CHOICE™ Adapters for single, row or column pipetting

Format	Net Volume	Quality	Box / PU	Order Number
96	1000 µl	PCR-certified	16	OL3811-25-539-N
		Pre-sterilized and PCR-certified		OL3811-25-639-S
		With filter (PCR-certified, pre-sterilized)		OL3811-25-939-F
	250 µl	PCR-certified	10	OL3811-25-537-N
		Pre-sterilized and PCR-certified		OL3811-25-637-S
		(160 µl) With filter (PCR-certified, pre-sterilized)		OL3811-25-937-F
	50 µl	PCR-certified	10	OL3811-25-535-N
		Pre-sterilized and PCR-certified		OL3811-25-635-S
		(40 µl) With filter (PCR-certified, pre-sterilized)		OL3811-25-935-F
384*	50 µl	Standard	10	OL0008-25-575-N

CyBi®-FeliX Head R 96/250 µl

using CyBi®-RoboTipTrays for 96-channel pipetting

Format	Net Volume	Quality	Trays/ PU	Order Number
96 Deep Well	250 µl	Standard	48	OL3810-26-661
		Pre-sterilized and PCR-certified		OL3810-26-664
		With filter (PCR-certified, pre-sterilized)		OL3810-26-669

using liquid handling adapters for single, row or column pipetting

Format	Net Volume	Quality	Box/ PU	Order Number
96	250 µl	PCR-certified	10	OL3811-25-537-N
		Pre-sterilized and PCR-certified		OL3811-25-637-S
		(160 µl) With filter (PCR-certified, pre-sterilized)		OL3811-25-937-F
	50 µl	PCR-certified	10	OL3811-25-535-N
		Pre-sterilized and PCR-certified		OL3811-25-635-S
		40 µl With filter (PCR-certified, pre-sterilized)		OL3811-25-935-F

CyBi®-FeliX Head R 96/60 µl

using CyBi®-RoboTipTrays for 96-channel pipetting

Format	Net Volume	Quality	Trays/ PU	Order Number
96	60 µl	Standard	48	OL3810-26-441
		Pre-sterilized and PCR-certified		OL3810-26-444
	40 µl	Standard		OL3810-26-631
		Pre-sterilized and PCR-certified		OL3810-26-634

using liquid handling adapters for single, row or column pipetting

Format	Net Volume	Quality	Box/ PU	Order Number
96	50 µl	PCR-certified	10	OL3811-25-535-N
		Pre-sterilized and PCR-certified		OL3811-25-635-S
	(40 µl)	With filter (PCR-certified, pre-sterilized)		OL3811-25-935-F

CyBi®-FeliX Head R 384/60 µl

using CyBi®-RoboTipTrays for 384-channel pipetting

Format	Net Volume	Quality	Trays/ PU	Order Number
384	60 µl	Standard	48	OL3810-26-241
		Pre-sterilized and PCR-certified		OL3810-26-244
	40 µl	Standard		OL3810-26-231

using Liquid Handling Adapters for single, row or column pipetting

Format	Net Volume	Quality	Box/ PU	Order Number
384*	50 µl	Standard	10	OL0008-25-575-N



CyBi®-Tip Box



CyBi®-RoboTipTray

* The CyBi®-TipBox for 384-well format (4.5 mm pitch) is racked with 192 tips for column or single pipetting.